


## Instruction for Use

# gastroplexVirus PLUS

## real time RT-PCR Kit

For qualitative *in-vitro* detection of RNA of Sapovirus and Astrovirus in clinical specimens, environmental and food samples.

<b>REF</b>	G01103-32	G01103-96
	32	96



gerbion GmbH & Co. KG  
Remsstr. 1  
70806 Kornwestheim  
Germany  
phone: +49 7154 806 20 0  
fax: + 49 7154 806 20 29  
e-mail: info@gerbion.com  
www.gerbion.com



## Index

1	Intended Use .....	3
2	Pathogen Information .....	3
3	Principle of the Test .....	3
4	Package Contents .....	4
5	Equipment and Reagents to be Supplied by User .....	4
6	Transport, Storage and Stability .....	5
7	Important Notes .....	5
8	General Precautions .....	5
9	Sample Material .....	5
10	Sample Preparation .....	5
11	Control RNA .....	6
12	Real time RT-PCR .....	7
12.1	Important Points Before Starting: .....	7
12.2	Procedure .....	7
12.3	Instrument Settings .....	9
13	Data Analysis .....	11
14	Assay Validation .....	13
15	Limitations of the method .....	13
16	Troubleshooting .....	13
17	Kit Performance .....	15
17.1	Diagnostic Sensitivity and Specificity .....	15
17.2	Analytical Sensitivity .....	15
17.3	Analytical Specificity .....	17
18	Abbreviations and Symbols .....	18
19	Literature .....	18

## 1 Intended Use

The gastroplexVirus PLUS real time RT-PCR Kit is an assay for the detection of RNA of *Sapovirus* and *Astrovirus* in clinical specimens (e.g. stool samples, vomit) real time PCR microplate systems.

## 2 Pathogen Information

Acute gastroenteritis is a worldwide major cause of morbidity and mortality. Enteric viruses are the major pathogens for gastroenteritis especially in children. Noro-, Rota-, Adeno-, Sapo- and Astroviruses are the most important viral pathogens.

**Astroviruses** are single stranded RNA (ssRNA) Viruses belonging to the family of Astroviridae. Diarrhoea is the most prevalent symptom of an Astrovirus-associated gastroenteritis, but also concomitant symptoms like vomiting and fever are described. In industrial countries, the incidence is 2-9%, especially in young children of under 2 years. Most relevant are the serotypes 1-5 of 8 serotypes known to date. The infection occurs by contaminated food and water or through the fecal-oral pathway.

**Sapoviruses** belong to the family of Calciviridae. Along with Noroviruses, Sapoviruses are the most common pathogens causing gastroenteritis worldwide. Although the highest incidence of Sapovirus infections is in young children under 5 years old, Sapovirus-associated gastroenteritis also occurs in adults. Clinical symptoms are similar to Norovirus infections like diarrhoea, vomiting, and fever, but the symptoms are milder.

To date less epidemiological studies are available and due to less sensitive diagnostic methods Sapoviruses were seldomly diagnosed.

## 3 Principle of the Test

The gastroplexVirus PLUS real time RT-PCR Kit contains specific primers and hydrolysis probes for the detection of the RNA of *Sapovirus* and *Astrovirus* in clinical specimens after the extraction of RNA from the sample material. The reverse transcription (RT) of viral RNA to cDNA and the subsequent amplification of virus specific fragments are performed in a one-step RT-PCR. The amplification can be detected when specific probes are hydrolysed by the Polymerase. The emitted fluorescence is measured in the FAM (*Sapovirus*) and ROX (*Astrovirus*) channel.

Furthermore, the gastroplexVirus PLUS real time RT-PCR Kit contains a Control RNA, which is detected in a second amplification system. Added during RNA

extraction, the Control RNA allows not only for the detection of RT-PCR inhibition but also detects possible mistakes during RNA extraction. This greatly reduces the risk of false-negative results. The Control RNA can also be used solely as Internal Control by adding it directly to the Mastermix. The fluorescence of the Control RNA is measured in the VIC®/HEX/JOE™/TET channel.

## 4 Package Contents

The reagents supplied are sufficient for 32 or 96 reactions respectively.

Table 1: Components of the gastroplexVirus PLUS real time RT-PCR Kit.

Label	Lid Colour	Content	
		32	96
Reaction Mix	yellow	1 x 506 µl	2 x 759 µl
Enzyme	blue	1 x 6.4 µl	1 x 19.2 µl
Positive Control <i>Sapovirus, Astrovirus</i>	red	1 x 50 µl	1 x 100 µl
Negative Control	green	1 x 50 µl	1 x 100 µl
Control RNA	colourless	1 x 160 µl	2 x 240 µl

## 5 Equipment and Reagents to be Supplied by User

- RNA isolation kit (e.g. NukEx Pure RNA/DNA, gerbion Cat. No. G05004 or NukEx Mag RNA/DNA, gerbion Cat. No. G05012)
- PCR grade Water
- Sterile microtubes
- Pipets (adjustable volume)
- Sterile pipet tips with filter
- Table centrifuge
- Vortexer
- Real time PCR instrument
- Optical PCR reaction tubes with lid
- Optional: Liquid handling system for automation
- Optional: VLP-RNA (gerbion Cat. No. G07008)

## 6 Transport, Storage and Stability

The gastroplexVirus PLUS real time RT-PCR Kit is shipped on dry ice or cool packs.

All components must be stored at maximum -18°C in the dark immediately after receipt. Do not use reagents after the date of expiry printed on the package. Up to 20 freeze and thaw cycles are possible.

For convenience, opened reagents can be stored at +2-8°C for up to 6 months. Protect kit components from direct sunlight during the complete test run.

## 7 Important Notes

- The gastroplexVirus PLUS real time RT-PCR must be performed by qualified personnel only.
- Good Laboratory Practice (GLP) has to be applied.
- Clinical samples must always be regarded as potentially infectious material and all equipment used has to be treated as potentially contaminated.

## 8 General Precautions

- Stick to the protocol described in the Instruction for Use.
- Set up different laboratory areas for the preparation of samples and for the set up of the RT-PCR in order to avoid contaminations.
- Pipettes, tubes and other materials must not circulate between those different laboratory areas.
- Always use filter tips.
- Regularly decontaminate equipment and benches with ethanol-free decontaminant.
- Do not combine gastroplexVirus PLUS real time RT-PCR Kit components of different lot numbers.

## 9 Sample Material

Starting material for the gastroplexVirus PLUS real time RT-PCR is the nucleic acid isolated from clinical specimens (e.g. stool samples, vomit), environmental or food samples.

## 10 Sample Preparation

The gastroplexVirus PLUS real time RT-PCR is suitable for the detection of *Sapovirus* and *Astrovirus* in clinical specimens (e.g. stool samples, vomit), environmental or food samples isolated with suitable isolation methods.

Commercial kits for RNA isolation such as the following are recommended:

- **NukEx Pure RNA/DNA**, gerbion Cat. No. G05004
- **NukEx Mag RNA/DNA**, gerbion Cat. No. G05012

**Important:** In addition to the samples always run a “water control” in your extraction. Treat this water control analogous to a sample.

Comparing the amplification of the Control RNA in the sample to the amplification of the internal control in the water control will give insights on possible inhibitions of the real time RT-PCR. Furthermore, possible contaminations during RNA extraction will be detectable.

**Please note the chapter ,Control RNA‘.**

If the real time RT-PCR is not performed immediately, store extracted RNA according to the instructions given by the RNA extraction kit’s manufacturer.

## 11 Control RNA

A Control RNA is supplied and can be used as extraction control or only as inhibition control. This allows the user to control the RNA isolation procedure and to check for possible real time RT-PCR inhibition.

The Virus-Like Particles (VLP-RNA) are not supplied, but must be added to the clinical or environmental samples directly. VLP-RNA can be used as patient-side extraction control and in automated extraction systems, when pipetting of the Control RNA to the first buffer of (e.g. binding buffer) of the respective extraction kit is not possible due to extraction instrument specifications.

### **RNA isolation from clinical specimens (e.g. stool samples, vomit), environmental and food samples**

#### a) Control RNA or VLP-RNA used as Extraction Control:

gastroplexVirus PLUS real time RT-PCR Control RNA or VLP-RNA is added to the RNA extraction.

Add 5 µl Control RNA or VLP-RNA per extraction (5 µl x (N+1)). Mix well. Perform the RNA isolation according to the manufacturer’s instructions. Please follow protocol A.

#### **The Control RNA must be added to the Lysis Buffer of the extraction kit.**

#### b) Control RNA used as Internal Control of the real time PCR:

If only inhibition will be checked please follow protocol B.

## 12 Real time RT-PCR

### 12.1 Important Points Before Starting:

- Please pay attention to the 'Important Notes' on page 5.
- Before setting up the real time RT-PCR familiarise yourself with the real time PCR instrument and read the user manual supplied with the instrument.
- The programming of the thermal profile should take place before the PCR set up.
- In every PCR run a Positive Control and one Negative Control should be included.
- Before each use, all reagents - except the Enzyme - should be thawed completely at room temperature, thoroughly mixed (do NOT vortex the Reaction Mix but mix by pipetting up and down repeatedly), and centrifuged very briefly.

### 12.2 Procedure

If the Control RNA or VLP-RNA is used to control both, the real time RT-PCR and the RNA isolation procedure, please follow protocol A. If the Control RNA is solely used to detect possible inhibition of the real time RT-PCR, please follow protocol B.

#### Protocol A

**The Control RNA or VLP-RNA was added during RNA extraction (see chapter 'Control RNA'). In this case, prepare the Master Mix according to Table 2.**

The Master Mix contains all of the components needed for RT-PCR except the sample. Prepare a volume of Master Mix for at least one sample more than required, in order to compensate for pipetting inaccuracy.

Table 2: Preparation of the Master Mix

Volume per Reaction	Volume Master Mix
15.8 µl Reaction Mix	15.8 µl x (N+1)
0.2 µl Enzyme	0.2 µl x (N+1)

**Protocol B**

The Control RNA is used for the control of the real time RT-PCR only (see chapter ,Control RNA'). In this case, prepare the Master Mix according to Table 3.

The Master Mix contains all of the components needed for real RT-PCR except the sample. Prepare a volume of Master Mix for at least one sample more than required, in order to compensate for pipetting inaccuracy.

**Important: Dilute the Control RNA 1:10 in PCR-grade dH<sub>2</sub>O (e.g. 1 µl Control RNA + 9 µl PCR grade Water before adding it to the Master Mix.**

Table 3: Preparation of the Master Mix

Volume per Reaction	Volume Master Mix
15.8 µl Reaction Mix	15.8 µl x (N+1)
0.2 µl Enzyme	0.2 µl x (N+1)
0.2 µl Control RNA* (diluted 1:10)	0.2 µl x (N+1)*

\*The increase in volume caused by adding the Control RNA is not taken into account when preparing the PCR assay. The sensitivity of the detection system is not impaired.

**Protocol A and B: real time RT-PCR set up**

- Place the number of optical PCR reaction tubes needed into the respective tray of the real time PCR instrument.
- Pipet **16 µl** of the Master Mix into each optical PCR reaction tube.
- Add **4 µl** of the eluates from the RNA isolation (including the eluate of the water control), the Positive Control, and the Negative Control to the corresponding optical PCR reaction tube (Table 4).
- Close the optical PCR reaction tubes immediately after filling in order to reduce the risk of contamination.



Table 4: Preparation of the real time PCR

Component	Volume
Master Mix	16.0 $\mu$ l
Sample	4.0 $\mu$ l
Total Volume	20.0 $\mu$ l

### 12.3 Instrument Settings

For the real time RT-PCR use the thermal profile shown in Table 5.

Table 5: real time RT-PCR thermal profile

Description	Time	Temperature	Number of Cycles
<i>Reverse Transcription</i>	10 min	45°C	1
<i>Initial Denaturation</i>	5 min	95°C	1
<i>Amplification of cDNA</i>			
Denaturation	10 sec	95°C	45
Annealing and Extension	40 sec	60°C	
	Aquisition at the end of this step		

Dependent on the real time instrument used, further instrument settings have to be adjusted according to Table 6.

Table 6: Overview of the instrument settings required for the gastroplexVirus PLUS real time RT-PCR Kit.

Real time PCR Instrument	Parameter	Detection Channel	Notes		
LightCycler 480II			Colour Compensation Kit Multiplex 1 (G070MP1-cc) required		
			Melt Factor	Quant Factor	Max Integration Time (sec)
	Sapovirus	FAM (465-510)	1	10	1
	Astrovirus	ROX (533-610)	1	10	2
	Control RNA	HEX (533-580)	1	10	2
	-	CY5 (618-660)	1	10	3
Stratagene Mx3000P / Mx3005P	Sapovirus	FAM	Gain 8		Reference Dye: None
	Astrovirus	ROX	Gain 1		
	Control RNA	HEX	Gain 1		
	-	Cy5	Gain 4		
ABI 7500	Sapovirus	FAM	Option Reference Dye ROX: NO		
	Astrovirus	ROX			
	Control RNA	JOE			
	-	Cy5			
Rotor-Gene Q, Rotor-Gene 3000 Rotor-Gene 6000	Sapovirus	Green	Gain 5		
	Astrovirus	Orange	Gain 5		
	Control RNA	Yellow	Gain 5		
	-	Red	Gain 5		

### 13 Data Analysis

The *Sapovirus* specific amplification is measured in the FAM channel and the *Astrovirus* specific amplification in the ROX channel.

The amplification of the Control RNA is measured in the VIC®/HEX/JOE™/TET channel.

Following results can occur:

- **A signal in the FAM or ROX channels is detected:  
The result is positive, the sample contains viral RNA.**  
In this case, detection of a signal of the Control RNA in the VIC®/HEX/JOE/TET channel is inessential, as high concentrations of cDNA may reduce or completely inhibit amplification of the Control RNA.
- **No signal in the FAM or ROX channels, but a signal in the VIC®/HEX/JOE/TET channel is detected:  
The result is negative, the sample does not contain viral RNA.**  
The signal of the Control RNA excludes the possibilities of RNA isolation failure (in case the Control RNA is being used as an Extraction Control) and/or real time RT-PCR inhibition. If the  $C_T$  value of a sample differs significantly from the  $C_T$  value of the water control, a partial inhibition occurred, which can lead to negative results in weak positive samples (see 'Troubleshooting').
- **Neither in the FAM and ROX channels nor in the VIC®/HEX/JOE/TET channel a signal is detected:  
A diagnostic statement cannot be made.**  
The RNA isolation was not successful or an inhibition of the RT-PCR has occurred. In case the Control RNA was added during RNA isolation and not directly to the PCR Master Mix, the Negative Control is negative in both channels.

**Figure 1** and **Figure 2** show examples for positive and negative real time RT-PCR results.




Figure 1: The positive sample shows virus-specific amplification in the FAM channel, whereas no fluorescence signal is detected in the negative sample.




Figure 2: The positive sample as well as the negative sample show a signal in the Control RNA specific VIC®/HEX/JOE/TET channel. The amplification signal of the Control RNA in the negative sample shows, that the missing signal in the virus specific FAM channel is not due to RT-PCR inhibition or failure of RNA isolation, but that the sample is a true negative.

## 14 Assay Validation

Set a threshold as follows:

### Negative Controls

All negative controls should be below the threshold. If there is a potential contamination (appearance of a curve in the negative control or a cluster of curves in specimens at high  $C_T$  – for example above 36), results obtained are not interpretable and the whole run (including extraction) has to be repeated.

### Positive Controls

All the positive controls must show a positive (i.e. exponential) amplification curve. The positive controls must fall below a  $C_T$  of 30.

### Internal Controls

All internal controls must show a positive (i.e. exponential) amplification curve. The internal control must fall below a  $C_T$  of 33. If the internal control is above  $C_T$  34, this points to a purification problem or a strong positive sample that can inhibit the IC. In the latter case, the assay is valid. If a water control run is performed, the IC must fall below a  $C_T$  of 33.

## 15 Limitations of the method

The results must always be considered in relation to the clinical symptoms. Therapeutical consequences should be made in consideration of clinical data. A negative test result does not exclude a *Sapovirus* or *Astrovirus* infection.

## 16 Troubleshooting

The following troubleshooting guide is included to help you with possible problems that may arise when performing a real time RT-PCR. If you have further questions, please do not hesitate to contact our scientists on [info@gerbion.com](mailto:info@gerbion.com).

### No fluorescence signal in the FAM or ROX channel of the Positive Controls

The selected channel for analysis does not comply with the protocol

Select the FAM channel for analysis of the Sapovirus specific amplification, the ROX channel for analysis of the Astrovirus specific amplification and the VIC®/HEX/JOE™/TET channel for the amplification of the Control RNA. Due to amplification both specific channels, amplification of the Internal Control can be inhibited in the Positive Control.

Incorrect configuration of the real time RT-PCR	Check your work steps and compare with 'Procedure' on page 7.
The programming of the thermal profile is incorrect	Compare the thermal profile with the protocol (Table 5, page 9).
Incorrect storage conditions for one or more kit components or kit expired	Check the storage conditions and the date of expiry printed on the kit label. If necessary, use a new kit and make sure kit components are stored as described in 'Transport, Storage and Stability', page 5.
<b>Weak or no signal of the Control RNA and simultaneous absence of a signal in the virus specific FAM channel or ROX channel.</b>	
real time RT-PCR conditions do not comply with the protocol	Check the real time RT-PCR conditions (page 6).
real time RT-PCR inhibited	Make sure that you use an appropriate isolation method (see 'Sample Preparation', page 5) and follow the manufacturer's instructions. Make sure that the ethanol-containing washing buffer of the isolation kit has been completely removed. An additional centrifugation step at high speed is recommended before elution of the RNA.
RNA loss during isolation process	In case the Control RNA was added before extraction, the lack of an amplification signal can indicate that the RNA isolation was not successful. Make sure that you use an appropriate isolation method (commercial kits are recommended) and stick to the manufacturer's protocol.
Incorrect storage conditions for one or more components or kit expired	Check the storage conditions and the date of expiry printed on the kit label. If necessary, use a new kit and make sure kit components are stored as described in 'Transport, Storage and Stability', page 5.
<b>Detection of a fluorescence signal in the FAM channel or ROX channel of the Negative Control</b>	
Contamination during preparation of the RT-PCR	Repeat the real time RT-PCR in replicates. If the result is negative in the repetition, the contamination occurred when the samples were pipetted into the optical PCR reaction tubes. Make sure to pipet the Positive Control last and close the optical PCR reaction tube immediately after adding the sample. If the same result occurs, one or more of the kit components might be contaminated. Make sure that work space and instruments are decontaminated regularly. Use a new kit and repeat the real time RT-PCR.

## 17 Kit Performance

### 17.1 Diagnostic Sensitivity and Specificity

During the validation study of the gastroplexVirus PLUS real time RT-PCR Kit 163 positive and 89 samples, negative for enteric pathogens were tested. The positive samples were also tested for the other pathogens to exclude unspecific reactions. The diagnostic sensitivity was found to be 100% and the diagnostic specificity 100%.

The positive predictive value was found to be 100%, the negative predictive value showed to be 100%.

Table 7: Overview of the amount of samples tested and the resulting sensitivity and specificity.

	gastroplexVirus PLUS positiv	gastroplexVirus PLUS negativ
Sapovirus positive	23	0
Sapovirus negative	0	229
Astrovirus positive	21	0
Astrovirus negative	0	231
<b>Sensitivity</b>	100%	
<b>Specificity</b>	100%	

### 17.2 Analytical Sensitivity

The limit of detection (LoD) of gastroplexVirus PLUS real time RT-PCR was determined using serial dilutions of in vitro transcripts (RNA Viruses) in nucleic acid stabilization buffer in a Stratagene Mx3000 real time PCR instrument. Total nucleic acids were extracted using NukEx Pure RNA/DNA according to the manufacturer's instructions. Each sample (200 µl of diluted nucleic acid) was supplemented with 5 µl Control-RNA prior to extraction. Total nucleic acids were eluted with 50 µl and 4 µl of the eluates were applied to the subsequent real time RT-PCR.

The LoD of gastroplexVirus PLUS real time RT-PCR for *Sapovirus and Astrovirus* is  $\geq 10$  copies per reaction each.

Table 8: Samples tested for the validation of the sensitivity of gastroplexVirus PLUS real time RT-PCR.

Sapovirus	Copies per Reaction	Expected Result	CT-value gastroplexVirus PLUS	Mean CT
Sapo 10-2	1.000.000	positive	23.34/ 22.71	23.03
Sapo 10-3	100.000	positive	26.93/ 26.58	26.76
Sapo 10-4	10.000	positive	30.57/ 30.27	30.42
Sapo 10-5	1.000	positive	33.85/ 34.38	34,12
Sapo 10-6	100	positive	38.11/ 37.64	37.88
Sapo 10-7	10	positive	40.36/ 40.79	40.58
Sapo 10-8	1	positive	no Ct/ no Ct	-

Astrovirus	Copies per Reaction	Expected Result	CT-value gastroplexVirus PLUS	Mean CT
Astro 10-2	1.000.000	positive	22.76/ 23.82	23.29
Astro 10-3	100.000	positive	25.75/ 25.56	25.66
Astro 10-4	10.000	positive	28.96/ 28.7	28.83
Astro 10-5	1.000	positive	31.48/ 31.81	31.65
Astro 10-6	100	positive	35.47/ 35.2	35.34
Astro 10-7	10	positive	39.81/ 38.77	39.29
Astro 10-8	1	positive	no Ct/ no Ct	-



### 17.3 Analytical Specificity

The specificity of gastroplexVirus PLUS real time RT-PCR was evaluated by *in silico* analysis and additionally by amplification of RNA and DNA of other relevant viruses and bacteria found in clinical samples.
















#### Results:

The gastroplexVirus PLUS real time RT-PCR showed a positive result for the samples containing *Sapovirus*, and *Astrovirus*, whereas samples containing other pathogens were reliably tested negative. The results are shown in Table 9.

Table 9: Bacterial and viral pathogens tested for the determination of the analytical specificity of gastroplexVirus PLUS real time RT-PCR.

Strain	Expected Result	Result
<i>Enterovirus 68</i>	negative	negative
<i>Coxsackievirus B3</i>	negative	negative
<i>Coxsackievirus A16</i>	negative	negative
<i>Coxsackievirus B5</i>	negative	negative
<i>Salmonella</i>	negative	negative
<i>Listeria monocytogenes</i>	negative	negative
<i>Escherichia coli</i>	negative	negative
<i>Campylobacter</i>	negative	negative
<i>Shigella</i>	negative	negative
<i>Yersinia</i>	negative	negative
<i>Norovirus G1</i>	negative	negative
<i>Norovirus GII</i>	negative	negative
<i>Rotavirus</i>	negative	negative
<i>Adenovirus</i>	negative	negative
<i>Sapovirus</i>	positive	positive
<i>Astrovirus</i>	positive	positive

## 18 Abbreviations and Symbols

cDNA	complementary Deoxyribonucleid Acid		Catalog number
RNA	Ribonucleid Acid		Contains sufficient for <n> test
PCR	Polymerase Chain Reaction		Upper limit of temperature
RT	Reverse Transcription		Manufacturer
	Reaction Mix		Use by YYYY-MM
	Enzyme		Batch code
	Positive Control		Content
	Negative Control		Consult instructions for use
	Control RNA		<i>In vitro</i> diagnostic medical device
			European Conformity

## 19 Literature

[1] Lothar Thomas, Labor und Diagnose: Indikation und Bewertung von Laborbefunden für die medizinische Diagnostik, 8. Auflage, 2012, TH-Books, ISBN-10: 3980521583