Gebrauchsinformation
Instruction Manual

NukEx Royal TS

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Index

Components	2
Abbreviations	2
Transport and Storage	2
Intended Use	2
Sample Material	2
Quality Control	3
Product Warranty	3
Introduction	3
Mode of Action	3
Equipment and Reagents to be Supplied by User	4
Important Notes	4
Procedure	5
Storage of Digestion Supernatants	6
Troubleshooting	6

Components

NukEx Royal TS

Nucleic Acid Release Reagent Tissue Shred

100 vials, each containing 300 µl NukEx Nucleic Acid Release Reagent and NukSand.

For 100 reactions.

Abbreviations

DNA	Desoxyribonucleic Acid
RNA	Ribonucleid Acid
PCR	Polymerase Chain Reaction
SNP	Single Nucleotide Polymorphis

Transport and Storage

The **NukEx Royal TS** is shipped on dry ice. It must be stored at -18 °C in the dark immediately after receipt. Do not use reagent after the date of expiry printed on the package. After initial usage, the reagent is stable for up to six months.

The reagent should not be thawed and frozen repeatedly.

Intended use

NukEx Royal TS can be used for the release of nucleic acid (DNA and RNA) from tissue samples (e.g. cartilage tissue, skin biopsies, adenoid tissue), ticks and insects through mechanical break down and subsequent enzymatic lysis. The supernatant obtained with **NukEx Royal TS** can be used for numerous molecular analyses.

G05001:

1.5 ml safe lock reaction tubes for manual use with a pestle (e.g. gerbion NukEx Pestle 1.5 ml Cat. No. G06006)

G05002:

2.0 ml screw cap tubes for automated use in e.g. FastPrep Cell Disrupter (qbiogene).

G05003:

2.0 ml safe lock reaction tubes for automated use in e.g. TissueLyser (Qiagen).

Sample material

Tissue samples like skin biopsies or adenoid tissue, ticks and insects.

Quality Control

In accordance with gerbion's ISO-certified Quality Management System, each lot of the **NukEx** Nucleic acid release reagent is tested against predetermined specifications to ensure consistent product quality.

Product Warranty

gerbion guarantees the performance of all products when used according to the instructions given in the Instruction Manual. The purchaser must determine the suitability of the product for its particular use. Should any product fail to perform satisfactorily due to any reason other than misuse, gerbion will replace it free of charge or refund the price. We reserve the right to change, alter, or modify any product to enhance its performance and design.

Introduction

For numerous molecular assays, the isolation of nucleic acid is required. For instance SNP analyses to proof the origin of animals, for the determination of animal species in meat analyses, or for the detection of viral and bacterial pathogens by PCR.

Mode of Action

NukEx Royal TS contains a special Sand which allows for the mechanical break down of tissue samples or insects and ticks. This break down can be done manually by using a pestle or automated in special instruments such as FastPrep Cell Disrupter (qbiogene) or TissueLyser (Qiagen). Subsequent enzymatic lyses of the samples with the NukEx Nucleic Acid Release Reagent (also contained in the NukEx Royal TS tubes) efficiently released nucleic acids. Crude supernatants can be employed directly in the subsequent molecular assay.

Depending on the assay to be performed, additional extraction of the nucleic acid with a commercial extraction kit (e.g. NukEx Pure RNA/DNA, gerbion Cat.No. G05004) is advisable.

Equipment and Reagents to be Supplied by User

General:

Laboratory equipment according to national safety instructions.

- Pipets (adjustable volume)
- Sterile pipet tips with filter
- Table centrifuge
- Block incubators or laboratory furnace
- dH₂O (PCR grade)

Sample preparation:

- Biopsy punch (app. 3-4 mm diameter; e.g. pfm AG, Cologne) or other instrument for sampling
- G05001: Pestle (e.g. gerbion NukEx Pestle 1.5 ml Cat. No. G06006)
- G05002/G05003: Instrument for homogenisation (e.g. FastPrep Cell disrupter or TissueLyser)

Optional:

- Liquid handling systems for automation
- Nucleic acid extraction kit (e.g. gerbion NukEx Pure RNA/DNA Cat. No. G05004)

Important Notes

- The **NukEx Royal TS** must be utilised by qualified personnel only.
- Good Laboratory Practice (GLP) has to be applied.
- Clinical samples must always be regarded as potentially infectious material and all equipment used has to be treated as potentially contaminated.

General Precautions

- Stick to the protocol described in the product insert.
- Always use filter tips.
- Regularly decontaminate equipment and benches with alcohol-free decontaminant.
- Use thermo-resistant labelling.

Procedure

1. Manual Tissue Break Down (G05001)

- Transfer sample into **NukEx Royal TS** tube
- Grind tissue with a pestle (e.g. NukEx Pestle 1.5 ml Cat. No. G06006)
- Ensure that samples are fully immersed in **NukEx Reagent**.
- Lock tube well.
- Incubate for 1 h at 60°C in a block incubator or laboratory furnace.
- Incubate sample for 10 min at 97 °C in a block incubator or laboratory furnace. This step is essential for the denaturation (inactivation) of the proteolytic enzymes. Unefficient inactivation can have negative impact on downstream applications!
- Let sample cool down for approx. 10 min at room temperature.
- Briefly spin down sample in a table centrifuge.
- Aliquots of the sample should be taken close to the surface in order to avoid the transfer of cell debris.

2. Automated Tissue Break Down (G05002 and G05003)

- Transfer sample into **NukEx Royal TS** tube
- Lock tube well.
- Place tube in homogenisation instrument (G05002: e.g. FastPrep Cell Disrupter; G05003: e.g. TissueLyser) and start program (the appropriate settings for the instruments are to be chosen by the user for each sample material used)
- Incubate for 1 h at 60°C in a block incubator or laboratory furnace.
- Incubate sample for 10 min at 97 °C in a block incubator or laboratory furnace. This step is essential for the denaturation (inactivation) of the proteolytic enzymes. Unefficient inactivation can have negative impact on downstream applications!
- Let sample cool down for approx. 10 min at room temperature.
- Briefly spin down sample in a table centrifuge.
- Aliquots of the sample should be taken close to the surface in order to avoid the transfer of cell debris.

Optional (G05001; G05002; G05003):

- Depending on the downstream application a purification of the supernatants is advisable. Usage of commercial extraction kits (e.g. gerbion NukEx Pure RNA/DNA Cat. No. G05004) is recommended.

Storage of Supernatants

Supernatants can be used immediately or can be stored under following conditions:

- up to 6 h at room temperature
- up to 24 h at +2 to +8 °C
- up to 6 months at -18 °C or lower temperature

In order to avoid cell debris being transferred, do not mix supernatants after thawing. Instead heat shortly at 60 °C and centrifuge if necessary.

Troubleshooting

The following troubleshooting guide is included to help you with possible problems that may arise in a subsequent PCR analysis. If you have further questions, please do not hesitate to contact our scientists on info@gerbion.com.

- **Neither the sample nor the Internal Control show a signal (inhibition):**
The heat inactivation of the proteolytic enzymes in **NukEx** was not effective. Heat the sample once more to 97 °C for 10 min. Repeat PCR analysis.
If necessary, purify nucleic acid from the supernatant by use of a commercial extraction kit (e.g. gerbion NukEx Pure RNA/DNA Cat. No. G05004) and repeat PCR analysis.
- **Weak signal of the Internal Control (partial inhibition):**
Components of the supernatant cause a partial inhibition of the PCR. Centrifuge the supernatant and dilute 1:3 in dilution buffer (e.g. **NukEx** Universal Dilution Buffer, G01014, gerbion). Repeat PCR with diluted supernatant.
- **The Internal Control still shows a weak signal after dilution:**
Purify nucleic acid from the supernatant by use of a commercial extraction kit (e.g. gerbion NukEx Pure RNA/DNA Cat. No. G05004) and repeat PCR analysis.